# PCR en temps réel - Light Cycller 480

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### Table des matières

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
The screenshot shows the 'Run Protocol' setup window. At the top, there are tabs for 'Run Protocol', 'Data', and 'Run Notes'. Below the tabs, there are fields for 'Detection Format' (SYBR Green I), 'Block Type' (96), 'Plate ID', and 'Reaction Volume' (20). The main section is divided into two parts: 'Programs' and 'Temperature Targets'.

**Programs Table:**

Program Name	Cycles	Analysis Mode
Program	1	None

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

ii. Subset

Create new analysis

Analysis Type \* Absolute Quantification

Subset \* All Samples

Program All Samples  
96er Format

Name

	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
A																								
B																								
C																								
D																								
E																								
F																								
G																								
H																								
I																								
J																								
K																								
L																								
M																								
N																								
O																								
P																								

✓ ✕

iii. Program

Create new analysis

Analysis Type \* Absolute Quantification

Subset \* All Samples

Program \* PCR

Name

	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
A																								
B																								
C																								
D																								
E																								
F																								
G																								
H																								
I																								
J																								
K																								
L																								
M																								
N																								
O																								
P																								

✓ ✕

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Create new analysis

Analysis Type \* Absolute Quantification

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

- **Monitoring PCR with the SYBR Green I Dye**

Generation of PCR products can be detected by measurement of the SYBR Green I fluorescence signal. SYBR Green I intercalates into the dsDNA helix. In solution, the unbound dye exhibits very little fluorescence; however, fluorescence (measured at 530 nm) is greatly enhanced (100-fold) upon binding to DNA due to conformational changes. Therefore, during PCR, the increase in SYBR Green I fluorescence is directly proportional to the amount of dsDNA generated. Since SYBR Green I dye is very stable it is the reagent of choice when measuring total DNA quantity.

The following are the basic steps of DNA detection by SYBR Green I during real-time PCR on the LightCycler® 480 System:

## 1. Programming

- Login to LightCycler® 480 Software
- Create a **New Experiment**
- The module **Experiment** is automatically opened. It consists of 3 folders:
  - **Run Protocol** for programming the experiment
  - **Data** for online data display
  - **Run Notes** optional

### Program - Setup

- Detection Format: **SYBR Green I**
- **Customize :** **Integration Mode : Dynamic**
- **Reaction Volume:** **20 µl**

### Program - Programs and Temperature Targets

Program Name	Cycles	Analysis Mode
Pre-incubation	1	None
Amplification	40	Quantification
Melting Curve	1	Melting Curves
Cooling	1	None

Pre-incubation				
Target (°C)	Acquisition Mode	Hold (hh:mm:ss)	Ramp Rate (°C/s)	Acquisitions (per °C)
95	None	00:05:00	4.4	

Amplification				
Target (°C)	Acquisition Mode	Hold (hh:mm:ss)	Ramp Rate (°C/s)	Acquisitions (per °C)
95	None	00:00:10	4.4	
55	None	00:00:15	2.2	
72	Single	00:00:08 (ex) Base #/25	4.4	
Melting Curve				
95	None	00 00:05		
65	None	00:01:00		
95	Continuous			10

Cooling				
Target (°C)	Acquisition Mode	Hold (hh:mm:ss)	Ramp Rate (°C/s)	Acquisitions (per °C)
40	None	00:01:00	1.5	

## 2. Subset Editor

Define a new subset that includes all wells that contain samples for the quantification assay. The MWP worksheet below is an example.

LightCycler 480 SYBR Green Master kit

	1	2	3	4	5	6	7	8	9	10	11	12	
A													A
B													B
C													C
D													D
E													E
F													F
G													G
H													H
	1	2	3	4	5	6	7	8	9	10	11	12	

	A 1 e 3
	B 2 e 3

	S6 1 e 6
	S5 1 e 5
	S4 1 e 4
	S2 1 e 3
	S1 1 e 2
	NTC

### 3. Sample Editor

The Sample Editor includes four folders (General, Abs Quant, Color Comp and Tm Calling) to edit general sample information and to define analysis specific sample information.

Sample information can be edited in these folders or by using the **Property Editor**.

- General: edit Sample Names and define Replicates
- Abs Quant: define Sample Types for Quantification and edit the respective Concentration for standard samples

#### Property Editor: For the definition of replicates

- Subset: User defined name
- Analysis Module: **General**
- Detection formats: SYBR Green I
- Sample Properties: Select individual subset groups and define them as **Replicates** of first well in row A (do not activate the original),

e.g. For Unknown A samples define as replicate of A1

**Apply:** Use this button before proceeding to next definition.

**Property Editor**

**Samples**

Subset: Quant Control Kit

	1	2	3	4	5	6	7	8	9	10	11	12
A												
B												
C												
D												
E												
F												
G												
H												

**Analysis**

Module: General

**Detection formats**

Filter combination: SYBR Green I

**Sample properties**

☐ Sample name

☒ Replicate: A1

☐ Sample ID

☐ Notes

OK Apply Cancel

### Property Editor: for definition of Standards and Concentrations

- Subset: User defined name
- Analysis Module: **Absolute Quantification**
- Detection formats: SYBR Green I
- Sample Properties: Activate the respective samples in the grid, define their **Sample Type** as **Standard** and edit the respective **Concentration**.
- Apply: Use this button and go on with defining the next standard dilution.

**Property Editor**

**Samples**

Subset: Quant Control Kit

	1	2	3	4	5	6	7	8	9	10	11	12
A												
B												
C												
D												
E												
F												
G												
H												

**Analysis**

Module: Absolute Quantification

**Detection formats**

Filter combination

SYBR Green I

**Sample properties**

☐ Target name

☒ Sample type: Standard

☒ Concentration: 100

OK Apply Cancel

#### 4. Pipetting

- Prepare PCR mixes and pipet into the MWP using template in part 3.

Component	1 Reaction	µl for <u>10</u> Reactions
Water, PCR grade (vial 2, colorless cap)	9,8 µl	98 µl
PCR primer mix (100 nM)	0,2 µl	2 µl
Master Mix, 2 x conc. (vial 1, green cap)	10 µl	100 µl
<b>TOTAL VOLUME OF MIX</b>	<b>18 µl</b>	<b>18 µl</b>
Template Volume	2 µl	2 µl
<b>TOTAL REACTION VOLUME</b>	<b>20 µl</b>	<b>20 µl</b>

- Cover the MWP with a sealing foil; use an adhesive seal applicator for correct sticking.
- Centrifuge the MWP at 1500 g for 2 minutes.

#### 5. Start Run

#### 6. Analysis

- **Absolute Quantification**

- **Analysis Type:** Absolute Quantification
- **Subset:** Subset from Part 3.
- **Program:** PCR

- When Analyses Window appears, select  . Standard curve will be plotted.

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Detection Format Profile: \_\_\_\_\_

Volume: \_\_\_\_\_

Program				
Program Names	Cycles		Analysis Mode	
Pre-incubation	1		None	
Amplification			Quantification	
Melting Curve	1		Melting Curves	
Cooling	1		None	
Temperature Targets				
Target (°C)	Acquisition Mode	Hold (hh:mm:ss)	Ramp Rate (°C/s)	Acquisitions (per °C)
Pre-incubation				
95	None	00:05:00 *	4.4	-
Amplification				
	None		4.4	-
	None		2.2	-
	Single		4.4	-
Melting Curve				
	None			-
	None			-
	Continuous		-	5-10
Cooling				
40	None	00:00:10	1.5	

\* If high polymerase activity is required, the pre-incubation program can be extended to 10 minutes.

Component	1 Reaction	µl for ____ Reactions
Water, PCR grade (vial 2, colorless cap)	µl	
PCR primer mix	µl	
Master Mix, 2 x conc. (vial 1, green cap)	µl	
TOTAL VOLUME OF MIX	µl	
Template Volume	µl	
TOTAL REACTION VOLUME	µl	

# LightCycler® 480 Worksheet – 96 well block

## - Probes Master with Hydrolysis Probes-

File Name: \_\_\_\_\_

Date: \_\_\_\_\_

User: \_\_\_\_\_

Detection Format Profile: \_\_\_\_\_

Volume: \_\_\_\_\_

Excitation Channels used: ☐ 450 ☐ 483 ☐ 523 ☐ 558 ☐ 615Emission Channels used: ☐ 500 ☐ 533 ☐ 568 ☐ 610 ☐ 640 ☐ 670

Program				
Program Names	Cycles		Analysis Mode	
Pre-incubation	1		None	
Amplification			Quantification	
Cooling	1		None	
Temperature Targets				
Target (°C)	Acquisition Mode	Hold (hh:mm:ss)	Ramp Rate (°C/s)	Acquisitions (per °C)
Pre-incubation				
95	None	00:05:00 *	4.4	-
Amplification				
	None		4.4	-
	None		2.2	-
	Single		4.4	-
Cooling				
40	None	00:00:10	1.5	

\* If high polymerase activity is required, the pre-incubation program can be extended to 10 minutes.

Component	1 Reaction	µl for ____ Reactions
Water, PCR grade (vial 2, colorless cap)	µl	
PCR primer mix	µl	
Master Mix, 2 x conc. (vial 1, red cap)	µl	
TOTAL VOLUME OF MIX	µl	
Template Volume	µl	
TOTAL REACTION VOLUME	µl	

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