


# Operating Instructions

**LKB**  
WALLAC



# CompuGamma Gamma Counter

Models 1282-001 and -002

The LKB-Wallac 1282 CompuGamma is an automatic micro-computer controlled universal gamma counter which can count samples of gamma emitting isotopes with energies in the range 10 to 2000 KeV. The standard model functions in two modes, Mode 0 allows a spectrum plot of CPM values, Mode 1 allows measurement in two channels of CPM values for single and double labelled samples.

CompuGamma is operated by typing instructions with a teleprinter keyboard. Results are printed out using the teleprinter and instantaneous count information appears on the built-in display. Step by step operating instructions are given in this manual, including printout examples with explanations.

## CONTROL LETTERS

Control Letter	Operation	See Page
P	Allows parameter values to be set	4,6,8a
A	Starts automatic counting	5,7b,8b
O	Interrupts any operation and sets the instrument to the READY state	5b,7b,8b
L	Lists the stored parameter groups	3b
C	Allows a parameter group to be cleared from the memory	3b
M	Moves parameters from one group to a new one	3b
I	Lists the control letters	3a
D	Instructs the numeric display to show a particular quantity	12a
R	Allows counting to be resumed	10a
Q	Allows counting of one rack	11a
E	Interrupts counting of a sample in the measuring chamber and allows editing of parameters	10a

Note: Page 6b means page 6 right hand half, page 15a means page 15 left hand half

## FEATURE INDEX

Feature	Page
Coding and correction	6b,15a
Display selection	12a
Display - spectrum	12a
Editing - interrupt	10b
Error numbers	16,17
Error - typing correction	6b
Factor - multiplying by	8a
Half life correction	8a
Interrupt	10-11
Jump to new line	6b
Level no. v. energy	13
Loading sample	5a
Low count reject	11b
Master Clear	15
Moving parameters	3b
Quick assay	11a
Peak integration	5b
Plug - code,pass,stop	14
Power failure	15b
Printout selection	9
Tabulation in printout	8a,9a
Text in printout	8a,9a
Window setting	14,15

Overprint indicates important points, and in the figures information typed by the user

# INTRODUCTION

## USING THE MANUAL

This page shows the use of two of the control letters I and L. The different operating modes are described with examples in the pages following: Mode 0, Spectrum Plot, on pages 4 and 5; Mode 1, CPM counting, on pages 6 and 7 with single label samples, and on page 8 with double label samples, decay correction etc. Other features are described in the rest of the manual from page 9 onwards. If Mode 0 is required start reading from page 4, if Mode 1 start from page 6.

## SWITCHING ON

If the instrument and output device are not in operation switch them on, and if a teleprinter is being used, set it for on-line operation.

## CONTROL LETTER LISTING

The control letters can be listed by typing I and RETURN as shown in the example

```
READY ->? I
          - Type I and press RETURN to
          list control letters

PARAMETER EDIT  C) CLEAR GROUP  L) LIST GROUPS  D) DISPLAY
AUTOMATIC OP    Q) QUICK ASSAY  E) EDIT OFF     R) RESUME
```

## EDITING OF STORED PARAMETER GROUPS

The number of parameter groups that are stored and the amount of memory available can be printed out by typing L. A group can be deleted from the memory by typing C. Parameters in one group can be moved to a new group by using control letter M; they cannot be moved to an existing group without that group being cleared.

```
READY -? L
          - Type L and press RETURN
          Stored groups listed
GR MODE SIZE% ID
1 1 2.6 BG-TEST
2 1 2.6 CS-137 EFF
3 0 1.6 SPECTRUM
20 1 3.3 DOUBLE

FREE MEMORY 89.9 %
          - Percentage of the memory free

READY -? C
PARAMETER GROUP (1-99) -? 1 - Type C to start clearing a group
CLEAR: -? 1 - Give group number first time
PARAMETER GROUP (1-99) -? 1 - Give number second time. If this
is a different number than the
first ERROR 75 will be printed
The group is cancelled and the
free memory increases

FREE MEMORY 92.5 %

READY -? M
MOVE FROM - Type M to begin move operation
PARAMETER GROUP (1-99) -? 2 - Give the number of the group
TO -? 4 - Give a new group number. If a
PARAMETER GROUP (1-99) -? 4 - group already exists with this
number the parameters will not
be moved
FREE MEMORY 89.9% - Memory available decreases as
new group is stored
Note: End every line by pressing
RETURN
```

# MODE 0 SPECTRUM PLOT

## STANDARD PARAMETERS FOR MODE 0

CompuGamma includes standard values which can be set for a parameter group, as shown in the example. In addition the allowed range of values for each parameter is also given. To begin parameter setting type P and press RETURN.

READY->? P			
PARAMETER GROUP	(1-99)	->? 10	- Type P and press RETURN
CREATE GROUP:	10		- Give a group number (1-99)
			- Shows a new group has been created
ID(MAX 20 CHAR)		->? SPECTRUM	- Give group name (Max.20 characters)
MODE(0-1)		->? 0	- Select MODE 0
LINE(1-13)			- Press RETURN list standard param.
1 LISTING	N		- Y-Yes,parameter listing required
2 TIME	10		- N-No listing required
3 COUNTS	900000		- Count time(1-99 999) secs
4 BACKGROUND	0		- Maximum no.of counts (0-9 999 000)
5 PEAKPRESET	0		- Preset value subtracted from each counting window (0-999 999)CPM
6 START CH	0		- Base value for peaks (0-99 999)CPM
7 WINDOW WIDTH	1		- Starting level no. (1-255)
8 STEP	1		- No.of levels per window (1-256)
9 STOP CH	256		- No.of levels from first level in one channel to first in next (1-255)
10 RAW DATA	N		- Level no.at which counting ends
11 REPEAT	1		- Y-Yes,results printed before plot
12 TTY LINE	72		- N-No,only plot required
13 FIRST POS	1		- No.of times sample counted (1-999)
LINE(1-13)		->? /	- Max.width of printout(72-132)
READY->?			- Sequence no. of first sample (0-999)
			- Jump to the READY state


## SPECTRUM PLOT PARAMETER SETTING

Parameter values can be changed and the new parameters stored with their own group number. The number of groups that can be stored depends on their sizes: at least 20 Mode 0 groups can be stored; Mode 1 groups are larger so less of them can be stored.

READY->? P			
PARAMETER GROUP	(1-99)	->? 10	- Type P to set parameters
EDIT GROUP:	10		- Select the parameter group number
ID(MAX 20 CHAR)	SPECTRUM->?		
MODE(0-1)	0	->?	
LINE(1-13)		->? 5	- Select first line to be changed
5 PEAKPRESET	0	->? 2000	- Peak base level set at 2000CPM
6 START CH	0	->? 150	- Spectrum begins from level no.150
7 WINDOW WIDTH	1	->?	
8 STEP	1	->?	
9 STOP CH	256	->? 190	- Spectrum ends at level no.190
10 RAW DATA	N	->? Y	- A printout of results is required
11 REPEAT	1	->?	
12 TTY LINE	72	->? 120	- Max.width of plot 120 characters
13 FIRST POS	1	->?	
LINE(1-13)		->? /	- Jump to READY state
READY ->?			Note: Press RETURN to end each line.

## LOADING SAMPLES ON THE CONVEYOR

Load samples - make sure that the racks are loaded so the two slots are towards the user, see the figure below. Each set of samples should have its own group selection plug, see page 14. Put a STOP plug after the last sample to stop the instrument. All sample holders should have a colour code ring, this fits into the top of the holder.



## AUTOMATIC SPECTRUM PLOT

Type A and give the group number to start automatic counting. Results are first listed then the spectrum is plotted out as shown in the example in the next column.

## STOPPING AUTOMATIC OPERATION

Type O and press RETURN to interrupt automatic counting and cause a jump to the READY state or use a STOP code plug in the last rack to be counted.

```

READY ->? A
PARAMETER GROUP (1-99) ->? 10
ID: SPECTRUM
MODE 0
- Type A to begin automatic counting
- Select parameter group
- Identifier and mode are printed automatically
  
```

LEVEL NO	CH 1	CH 2	CH 3	CH 4	CH 5	TIME	ETIME
150-155	770	629	931	552	610	10	0.00
155-160	777	610	571	616	571	10	0.01
160-165	668	713	847	1111	1342	10	0.01
165-170	1060	1195	1394	1439	1695	10	0.01
170-175	2439	3185	6104	9518	15840	10	0.01
175-180	20863	19753	15846	11458	6095	10	0.02
180-185	3178	1310	571	315	250	10	0.02
185-190	141	180	122	122	160	10	0.02

SPECTRUM PLOT: POS 1  
In the spectrum plot the CPM at each level number is given and the percentage that this is of the total CPM

LEVEL NO	% OF TOTAL CPM	CPM 0	5200	10400	20800
150-151	0.32	770 I	10400	20800	
151-152	0.43	629 I	10400	20800	
152-153	0.63	931 I	10400	20800	
153-154	0.37	552 I	10400	20800	
154-155	0.41	610 I	10400	20800	
155-156	0.53	777 I	10400	20800	
156-157	0.41	610 I	10400	20800	
157-158	0.39	571 I	10400	20800	
158-159	0.42	616 I	10400	20800	
159-160	0.39	571 I	10400	20800	
160-161	0.45	668 I	10400	20800	
161-162	0.48	713 I	10400	20800	
162-163	0.57	847 I	10400	20800	
163-164	0.75	1111 I	10400	20800	
164-165	0.91	1342 I	10400	20800	
165-166	0.72	1060 I	10400	20800	
166-167	0.81	1195 I	10400	20800	
167-168	0.95	1394 I	10400	20800	
168-169	0.98	1439 I	10400	20800	
169-170	1.15	1695 I	10400	20800	
170-171	1.65	2439 I	10400	20800	
171-172	2.16	3183 I	10400	20800	
172-173	4.14	6104 I	10400	20800	
173-174	6.46	9518 I	10400	20800	
174-175	10.74	15840 I	10400	20800	
175-176	14.15	20863 I	10400	20800	
176-177	13.40	19753 I	10400	20800	
177-178	10.75	15846 I	10400	20800	
178-179	7.77	11458 I	10400	20800	
179-180	4.13	6095 I	10400	20800	
180-181	2.16	3178 I	10400	20800	
181-182	0.89	1310 I	10400	20800	
182-183	0.39	571 I	10400	20800	
183-184	0.21	315 I	10400	20800	
184-185	0.17	250 I	10400	20800	
185-186	0.10	141 I	10400	20800	
186-187	0.12	180 I	10400	20800	
187-188	0.08	122 I	10400	20800	
188-189	0.08	122 I	10400	20800	
189-190	0.11	160 I	10400	20800	
150-190	100.00	147424	10400	20800	

Peak indicated by values with double lines (=) Integrated area under peak includes all values marked in colour here

PEAK 1  
170-181 75.52 114277  
- First peak greater than PEAKPRESET value: width of the peak, percentage of the total CPM within the peak, and CPM within the peak are given  
- Total CPM within this counting window 150 - 190

# MODE 1 CPM SINGLE LABEL

## STANDARD PARAMETERS FOR MODE 1

To set the standard parameters type P and press RETURN  
 CompuGamma includes standard values which can be set  
 for a parameter group, as shown in the example. In  
 addition the allowed range of values for each parameter is  
 also given. To begin parameter setting type P and press  
 RETURN.

READY- ? P	PARAMETER GROUP (1-99)	->? 1	
CREATE GROUP: 1		->? EXAMPLE	- Give a parameter group number (1-99)
ID (MAX 20 CHAR)			- Shows a new group has been created.
			- Give group name (max .20 characters)
MODE (0-1)		->? 1	Note: Do not use , ' / in the ID
LINE (1-20)		->?	- Select MODE 1
			- Press RETURN to list standard param.
1 LISTING	N		- Y-Yes, parameter listing required
2 TIME	60		- N-No listing
3 COUNTS 1	900000		- Count time (1-99 999) secs
4 COUNTS 2	900000		- Max.no. counts in ch.1(0-9 999 000)
5 LCR COUNTS	0		- Same for channel 2
6 LCR TIME	0		- Low count reject count(0-999 999)
7 BACKGROUND 1	0		- LCR time(0-999)secs
8 BACKGROUND 2	0		- Ch.1 preset background(0-999 999) CPM
9 WINDOW 1	035-102		- Ch.2 preset background(0-999 999) CPM
10 WINDOW 2	000-000		- Ch.1 window width(0-256)
11 HALF LIFE 1	0		- Ch.2 window width(0-256)
			- Half life, ch.1(0.001-9 999.99)hours
			- 0-no correction See Instr.Man.Sec.4
12 HALF LIFE 2	0		- Half life corr. for channel 2
13 FACTOR 1	1		- For multiplying RATIO(0.0-999 999.9)
14 FACTOR 2	1		- For multiplying RATIO2(0.0-999 999.9)
15 SUMMATION	N		- Y-Yes, add ch.1 and ch.2 results
			- N-No, do not add results
16 REPLICATE	1		- Number of replicates(1-99)
17 REPEAT	1		- No. of times sample counted(1-999)
18 LABELS(1/2)	1		- 1 - single label samples
			- 2 - double label samples
19 CODING	POS CODE		- Sample position and code
20 PRINT	-1-2-3-4-5-6-7-8		- Printout selection
LINE(1-20)		->? /	- Type / RETURN to jump to READY,
READY->?			Parameters are stored automatically
			Note: Type / 0 RETURN to go to LINE,
			and / xx RETURN to go to line xx

## CPM MEASUREMENT PARAMETER SETTING EXAMPLE 1

Parameter values can be changed: first type P and select  
 the parameter group then give the new parameter values,  
 as shown in the example. Always press RETURN at the  
 end of every line.

READY- ? P	PARAMETER GROUP (1-99)	->? 1	
EDIT GROUP: 1			- Parameter setting selected
ID (MAX 20 CHAR) EXAMPLE1->?			- Parameter group selected
			- Param.group created previously
MODE (0-1)	1	->? 6.6 2	- RETURN pressed as line unchanged
LINE (1-20)			- RETURN pressed
			- Line 2 selected, (6 was an error)
			Note: If a typing error occurs
			go back to the error by pressing
			DELETE and correct it
2 TIME	60	->? 40	- New count time set
3 COUNTS 1	900000	->? /16	- Jump to line 16
16 REPLICATE	1	->? 2	- Duplicates selected
17 REPEAT	1	->? 2	- Count to be repeated
18 LABELS(1/2)	1	->?	- RETURN pressed
19 CODING	POS-CODE		- Position and code for control
	000-	->? 1-BGND	samples and unknowns. No. of
	000-	->? 3-REFA	replicates for controls deter-
	000-	->? 5-UNKS	mined by coding not line 16
	000-	->?	- RETURN pressed to end coding
	-1-2-3-4-5-6-7-8		- New printouts selected
	->? 1-2-3-4-5-6-7-8-9		- RETURN pressed as no more coding
	->?		
20 PRINT			

### AUTOMATIC SAMPLE COUNTING

Type A and give the group number to start automatic counting.

POS CODE	ETIME	CTIME	DTIME%	COUNTS1	CPM1	ERR1%	RATIO1
001 BGN	0.04	40	0.6	165	249	7.8	
001 MEAN	0.05	40	0.6	166	251	7.8	
002	0.07	40	0.6	146	220	8.3	
002 MEAN	0.08	40	0.6	145	216	8.4	
003 REFA	0.09	40	1.4	18366	27697	0.7	
003 MEAN	0.11	40	1.4	18520	27933	0.7	
004	0.12	40	1.4	18326	27635	0.7	
004 MEAN	0.13	40	1.4	18072	27246	0.8	
005 UNKS	0.15	40	1.8	16422	24720	0.8	0.8947
005 MEAN	0.16	40	1.8	16312	24552	0.8	0.8887
006	0.17	40	1.8	16391	24672	0.8	0.8930
006 MEAN	0.19	40	1.8	16144	24294	0.8	0.8793
	0.20	40		24483	24483	0.6	0.8861
MEAN				24559		0.4	0.8889

### STOPPING AUTOMATIC OPERATION

Type O and press RETURN to interrupt automatic counting and cause a jump to the READY state or use a STOP code plug in the last rack to be counted.

LINE(1-20)	N	POS-CODE	RETURN
1 LISTING	40	001-BGND	1
2 TIME	900000	003-REFA	1
3 COUNTS 1	900000	005-UNKS	1
4 COUNTS 2	0	000-	1
5 LCR COUNTS	0	-1-2-3-4-5-6-7-8-9	1
6 LCR TIME	0		
7 BACKGROUND 1	0		
8 BACKGROUND 2	0		
9 WINDOW 1	035-102		
10 WINDOW 2	000-000		
11 HALF LIFE 1	0		
12 HALF LIFE 2	0		
13 FACTOR 1	1		
14 FACTOR 2	1		
15 SUMMATION	N		
16 REPLICATE	2		
17 REPEAT	2		
18 LABELS (1/2)	1		
19 CODING	1		
20 PRINT			

Note: RETURN must be pressed to end every line

### LOADING SAMPLES ON THE CONVEYOR

Load samples - make sure that the racks are loaded so the two slots are towards the user. Each set of samples should have its own group selection plug, see page 14. Put a STOP plug after the last sample to stop the instrument.

# MODE 0 SPECTRUM PLOT

## STANDARD PARAMETERS FOR MODE 0

CompuGamma includes standard values which can be set for a parameter group, as shown in the example. In addition the allowed range of values for each parameter is also given. To begin parameter setting type P and press RETURN.

LINE(1-13)	PARAMETER GROUP (1-99)	VALUE	DESCRIPTION
1 LISTING	N		- Y-Yes, parameter listing required
2 TIME	10		- N-No listing required
3 COUNTS	900000		- Count time(1-99 999) secs
4 BACKGROUND	0		- Maximum no.of counts (0-9 999 000)
5 PEAKPRESET	0		- Preset value subtracted from each counting window (0-999 999)CPM
6 START CH	0		- Base value for peaks (0-99 999)CPM
7 WINDOW WIDTH	1		- Starting level no. (1-255)
8 STEP	1		- No.of levels per window (1-256)
9 STOP CH	256		- No.of levels from first level in one channel to first in next (1-255)
10 RAW DATA	N		- Level no.at which counting ends
11 REPEAT	1		- Y-Yes, results printed before plot
12 TTY LINE	72		- N-No, only plot required
13 FIRST POS	1		- No.of times sample counted (1-999)
LINE(1-13)			- Max.width of printout(72-132)
READY->?			- Sequence no. of first sample (0-999)

## SPECTRUM PLOT PARAMETER SETTING

Parameter values can be changed and the new parameters stored with their own group number. The number of groups that can be stored depends on their sizes: at least 20 Mode 0 groups can be stored; Mode 1 groups are larger so less of them can be stored.

LINE(1-13)	PARAMETER GROUP (1-99)	VALUE	DESCRIPTION
5 PEAKPRESET	0	2000	- Peak base level set at 2000CPM
6 START CH	0	150	- Spectrum begins from level no.150
7 WINDOW WIDTH	1		
8 STEP	1		
9 STOP CH	256	190	- Spectrum ends at level no.190
10 RAW DATA	N	Y	- A printout of results is required
11 REPEAT	1		
12 TTY LINE	72	120	- Max.width of plot 120 characters
13 FIRST POS	1		
LINE(1-13)			- Jump to READY state
READY ->?			Note: Press RETURN to end each line.



# GENERAL

## PRINTOUT SELECTION

Information to be printed out can be selected according to individual requirements, it can be printed out in any order, special characters may be introduced, information can be repeated and can extend to as many printout lines as required. Each item is selected according to its reference number as shown in the accompanying table, for convenience the standard program includes a preset printout format. The following points should be noted:

1 Tabulation is possible by using the symbol >xx where xx is the number of characters from the left margin at which printing is to start, e.g. -5->20-6-means that printout column 6 starts 20 character positions in from the left margin. This assumes that printout column 5, and those before it, have ended before position 20 has been reached, if this is not so, the tabulation will be ignored in the printout. The tabulation instruction must be separated from the other printout options by hyphens (-) in the usual way.

2 Text can be included in the printout by putting it between apostrophes in print selection e.g. -5-'channel 1'-6-

3 The 's' in the printout column width listing means that an extra character place is reserved in the printout, this is either a space if the number is positive or a minus sign (-) if the number is negative. The 'SP' in the listing means that one extra space is included after every column. The numbers given are the total width of the printout column for that printout option.

4 If the number measured or calculated is too big for the column width assigned to it, the column will be filled with asterisks e.g. \*\*\*\*\*.

Ref No	Heading	Column width
0	Return and line feed	
1	Sequence number	xxxx SP 4
2	Code	xxxx SP 5
3	Elapsed time in hours	xxxx.xx SP 9
4	Counting time in seconds	xxxxxx SP 7
5	Dead time as a percentage of count time	xxxx.x SP 7
6	Counts in channel 1	COUNTS1 sxxxxxxx SP 9
7	CPM in channel 1	CPM1 sxxxxxxx SP 9
8	Error in ch1 CPM	ERR1% sxxx.x SP 7
9	Unk CPM ch1	RATIO1 sxx.xxxx SP 9
	Ref A mean CPM ch1 or for ref.sample A	
	Mean CPM ch2	
10	Ratio 1 x factor 1	FACTOR1 sxxxx.xx SP 9
11	Counts in ch2	COUNTS2 sxxxxxxx SP 9
12	CPM in ch2	CPM2 sxxxxxxx SP 9
13	Error in ch2 CPM	ERR2% sxxx.x SP 7
14	Unk CPM ch2	RATIO2 sxx.xxxx SP 9
	Ref B mean CPM ch2 or for ref.sample B	
	Mean CPM ch1	
	Mean CPM ch2	
15	Ratio 2 x factor 2	FACTOR2 sxxxx.xx SP 9
16	Counts 1 + counts 2	COUNTS# sxxxxxxx SP 9
17	CPM1 + CPM2	CPM# sxxxxxxx SP 9
18	Error in(ch1 + ch2) CPM	ERR# sxxx.x SP 7
19	Unk CPM1 + Unk CPM 2	RATIO# sxx.xxxx SP 9
	Ref A mean CPM1 + CPM2	
20	Ratio # x factor 1	FACTOR# sxxxx.xx SP 9

```

READY ->? A
PARAMETER GROUP (1-99) ->? 1
ID:TEST
MODE 1
- Counting is started

POS CODE    CTIME COUNTS1    CPM1    ERR1X
001 BGNB    60    21    21.8
002          60    28    18.9
          MEAN    25    14.3
          OObserve the count time
          and number of printout
          columns

003 REFA    60    27822    28182    0.6
004          60    28055    28421    0.6
          MEAN    28301    0.4
          OObserve the count time
          and number of printout
          columns

005 UNKS    60    23780    24056    0.6
          OObserve the count time
          and number of printout
          columns

*COUNTING INTERRUPTED*
PARAMETER GROUP 1
ID:TEST
MODE 1
POS 6
- Counting interrupted by
  typing letter E

EDIT ->? L
GR MODE SIZEX ID
2 1 2.2 DEMO
1 1 2.8 TEST
FREE MEMORY 95.0%
- After EDIT is printed other
  operations can be performed,
  e.g. Listing

EDIT ->? B
PARAMETER GROUP (1-99) ->? 1
EDIT GROUP:1
ID(MAX 20 CHAR) TEST ->?
MODE(O-1) 1 ->?
LINE(1-20) ->? 2
2 TIME 60 ->? 30
3 COUNTS 1 900000 ->? /20
20 PRINT -1-2-4-6-7-8
-1-2-4-6-7-8 ->? 1-2-4-6-7-8-9-10
->? /
- A new count time is set
- Printout options are changed

EDIT ->? R
*COUNTING RESUMED*
PARAMETER GROUP 1
ID:TEST
MODE 1
POS 6
- R is typed to resume counting
  with the same parameter group
  as when counting was
  interrupted (in this case with
  changed parameters)

POS CODE    CTIME COUNTS1    CPM1    ERR1X RATIO1    FACTOR    Observe
006          30    12014    24309    0.9 0.8589    0.86 count time
007          30    9107    18391    1.0 0.6498    0.65 and printout
008          30    9170    18518    1.0 0.6543    0.65 columns
009          30    6490    13080    1.2 0.4622    0.46
010          30    6565    13231    1.2 0.4675    0.47

```

## INTERRUPTING THE STANDARD CURVE

The control letter O can be typed at any time but if this happens during curve calculation the curve will be lost and no curve will be in the memory if an attempt is made to resume the assay, only CPM values will be calculated and printed.

Control letter E can be given at anytime but it will not interrupt calculation of the standard curve. The program only jumps to the EDIT state after the curve has been plotted.

## INTERRUPT FOR EDITING

Counting of a sample in the measuring chamber can be interrupted by typing E followed by pressing RETURN. Any of the control letters C)lear L)ist P)arameter R)esume and O)ff can be typed. When the editing procedure has been completed type R followed by RETURN; counting will resume in the normal way and will include any changes made during editing; it will begin with the sample in the measuring chamber which will be recounted.



## DISPLAY

Numerical information on the display is selected by typing D followed by a number, or a letter and a number, then pressing RETURN. Display possibilities are listed below.

Table of display information

D1	Position	
D2	Time	
D3	Elapsed time	
DWx	Window setting for ch x	where x is 1 - 5
DGx	Counts in ch x	
DCx	CPM in ch x	
D9	Scans each value and displays it for 4 secs	

## SPECTRUM DISPLAY

Part of the display also provides an indication of the energy distribution of gamma emission from the sample. Thirty-two light emitting diodes are arranged to cover the 256 energy levels, (eight levels for each diode,) they indicate by their intensity the spectral distribution of radiation from the sample. Each diode covers a segment of the energy range, the maximum energy for each is indicated by the energy level number marked for the individual diodes. Together they provide a rapid means of checking that the count window is set correctly for the sample.

## SETTING THE COUNTING WINDOW

The energy range of radiation detectable by the Computer is 10 - 2000 KeV; this range is divided logarithmically into 256 levels, each of which is given a number. The counting window is defined by giving the upper and lower level numbers. Once this window has been set only those particles with energies falling within in the range of the window will be counted by the microcomputer. There are various ways of deciding where to set the window limits:

### Preset window

There is a standard window setting given in Mode 1, this is from level number 35 to 102 for WINDOW 1 which is a good starting window for counting isotopes  $^{125}\text{I}$  and  $^{129}\text{I}$ , however the user may change this window and also WINDOW 2 which is set at 000-000 in the standard parameters.

### Spectrum Display


The 32 LEDs on the display show the peak position and allow approximate limits to be set.

### Spectrum Plot Mode 0

The main method for deciding which levels to set is by running a spectrum plot on the sample using Mode 0. This allows the optimum window to be set for a particular sample.


### LEVEL NUMBER/ENERGY PLOT

The figure below shows a plot of the log of energy against level number, from this the level number corresponding to a particular energy, or the energy range corresponding to a particular LFD may be obtained.



### ISOTOPE/LEVEL NUMBER PLOT

The figure below shows the peak positions in terms of the level number of some selected isotopes.



- Strongest peaks
- ▨ Intermediate peaks
- Weakest peaks

## CODE PLUGS

Automatic counting of racks can be controlled by the use of code plugs. A code plug occupies one position in a sample rack and when it reaches the counting position it gives information to the instrument. There are three main types of information:

### Normal numeric code plug

This plug has a number in the range 0-9. When the instrument detects one or two of these plugs it starts counting the samples that follow the plug(s) using the parameters in the group given by the number of the code plugs. Using two plugs, groups in the range 1-99 can be selected automatically. If LISTING has been set to Y parameters will be listed before results are printed. If N only the group number, identifier and mode will be printed out before count results.

Note: Not more than two numeric code plugs may be put in successive rack positions.


### Pass plug

When a pass plug is detected the samples following are not counted until another pass plug is detected at which point counting resumes normally. The message \* PASS will be printed out. Letter O should not be typed during the time between the detection of the first and last pass plugs. A pass plug must not be put after a numeric code plug with no samples in between.

### Stop plug

Counting can be stopped by inserting this plug. The rack will be driven out of the counting position and the instrument will jump to the READY state. The message \* END OF ASSAY will be printed out. A STOP plug will be ignored if it is in between two pass plugs.

Note: Any plugs in a rack the wrong way round will be ignored but the sample position number will be increased by ten after the rack has passed through the measuring position.



## CODING

The instrument recognises four special codes: BGND REFA REFB UNKS. No code is recognised after the UNKS code has been given, it will be treated as a name for unknown samples. If an unrecognised code is given before UNKS it will be treated as an unknown sample as far as possible e.g. if it comes between REFA and REFB it cannot be treated as a double label sample because the second reference sample has not been measured. It will be given the name typed in coding. If the next code is a special one e.g. REFB it will be treated in the normal way. The first three special codes can be given in any order but the background sample will only be subtracted from those samples coded after it. A second background sample can be coded before UNKS and will be subtracted from the unknown samples. The position numbers given to the samples should be in increasing order. Replicate samples must not be coded.

Editing of coding can be done in the following way:

If a line of coding is to be left out type - .

If only position number is to be changed give the new number only.

If only code is to be changed give - followed by the new code

If both number and code are to be changed then type both with the - in between in the normal way.

Note: To correct an error 74 be sure that the new position number given is greater than the position number before it in the coding and less than the one after it otherwise another error 74 will be indicated.

## POWER FAILURE

The CompuGamma is fitted with a standby power supply which protects it for up to 45 mins. against the loss of parameters if a power failure occurs. If power is off for a longer time all the parameters will be lost and the instrument will be in the MASTER CLEAR state when power returns. When power returns the instrument enters a power failure loop. This is indicated on the display by the numerical display showing a row of 8s and the other LEDs all being lit up except for the TIME indicator LED which goes on and off. This state will continue for about 15 minutes if uninterrupted, then the program will jump to the READY state unless it was in the automatic counting mode when power failure occurred in which case it will resume counting. This state can be interrupted in the following ways:

Typing C causes a jump to MASTER CLEAR and all the parameters will be lost.

Typing / causes a jump to the READY state but saves the parameters

Pressing RETURN causes the words POWER FAILURE to be printed then the program jumps to the READY state or continues counting if it was in the automatic counting mode when power was lost.

Note: The printer must always be reset by pressing the TERM READY button although results can still be printed out by the instrument if this has not been done.

## ERROR NUMBERS

If an error occurs during operation, the CompuGamma will inform the user by means of an error number. The meaning of each error number is described in the table following. Where necessary, instructions are given about the action that needs to be taken.

No.Error

- 6 The conveyor does not start. This may be caused by one or more motors not being at the starting point of their operational cycle, or by a sensor failure. Type letter O and press RETURN to drive the motors to their initial positions.
- 11 Operation cannot be resumed because there is a Stop plug, a Code plug or no sample in the first rack position to be counted.
- 41 Stabilisation failed during measurement.
- 42 Error when calculating CPM or counting error in CPM. The count time is very small or zero causing a division-by-zero error or an overflow error.
- 43 Error in spillover correction in double label assay. Ratios  $q(A)$  and  $q(B)$  are such that calculation of corrected CPM or its error leads to an error.
- 44 Error when calculating  $q(A)$ ,  $q(B)$ , their errors or  $1 - q(A) \cdot q(B)$  in a double label assay. Check the window settings, coding and reference samples to see if they are correct.
- 45 Overflow has occurred when summing repeat CPMs or their errors. This error can also result after errors 42, 43, and 47.
- 46 Overflow has occurred when summing CPM or CPM error values for the replicate mean value calculation. This can also result after errors 42, 43, 47 and 45.

No.Error

- 47 Error in half life correction. The ratio elapsed time/half life has become so large that calculation of corrected CPM or its error leads to overflow.
- 48 No double label nor half life correction has been made because one or both of the reference samples is missing.
- 49 Total sum of CPMs is zero (or nearly zero) leading to error in percent calculations (MODE 0).
- 50 Parameter group number given after A or Q command differs from the group number specified by the code plugs. Code plugs override the number given.
- 51 No group is stored with the group number specified by the code plugs. All samples up to the next code or stop plug are ignored.
- 53 Rack error. The rack has been put the wrong way round on the conveyor. The samples in the rack will not be counted but the sample positioning number will increase by ten for the next sample. Any plugs in the rack will be ignored.
- 57 Checksum error. An error has been introduced into the BASIC program in the non-permanent memory, the program is reloaded automatically.
- 58 Error discovered which prevents the instrument going to the READY state, instead it goes to the service state; call a service engineer.
- 59 Operation has stopped because the BASIC interpreter has detected an error while the program was running. Type the same control letter as was used to start the operation. If the same error occurs again (and also after Master Clear) it is an indication of a memory fault. This error message will be followed by another of the form ?xx Error in yyyy, then READY? will be printed out.



No.Error

- 61 Command is not recognised. Type I to list those control letters allowed.
- 62 No parameter group is stored with this number. Check list of parameter groups by typing L followed by pressing RETURN.
- 63 An integer number was expected as input and instead some other characters have been given.
- 64 Illegal characters given in the numeric data.
- 65 Numerical value given as input is outside the range allowed. See the Reference Card.
- 66 Instruction to jump to a parameter line cannot be followed because either the line number is not correct, or the mode has not been given so the number of lines is not defined.
- 67 An illegal window has been set: windows must be given in the format xxx-yyy where  $0 \leq xxx \leq yyy \leq 256$ .
- 68 Incorrect printout selection, illegal number given.
- 69 Incorrect printout selection, tabulation exceeds line width.
- 70 A number must follow the tab (>) character.
- 71 Illegal character in printout selection: text must be given inside apostrophes (').
- 72 Printout selection too long, more than 80 items.
- 73 Y(ES) or N(O) must be given.
- 74 Difference in position number in coding must be within the range 1 - 99. If the earlier position number is greater than the second then 1000 is added to the second and the test is made again. E.g. the difference between samples coded with position 997 and 003 is  $1003 - 997 = 6$ .

No.Error

- 75 Parameter group cannot be cleared because the group number given the second time (in response to the group number confirmation question) differs from that given first. Start clear procedure again from the READY state by typing C.
  - 80 Mode must be given when creating a new group.
  - 81 Memory size exceeded. Jump to the READY state and cancel some other parameter group from the memory.
  - 82 Cannot resume because no assay interrupted.
  - 83 Cannot resume because the mode in the parameter group has been changed.
  - 84 A checksum test has revealed an error in a parameter group in the memory. The numbers of the groups which have been lost are listed.
  - 91 One or more motors have been on too long (more than 9 minutes) so all motors have been stopped. Possible causes are:
    - a) An obstacle on the conveyor. Remove the obstacle and type letter O then press RETURN to drive the motor to the initial position.
    - b) The instrument has been left running in the A mode with no racks, or invalid racks, for more than 9 minutes. Continue as in a)
    - c) A motor is not working properly.
- Note: If the message RAM FAILURE xxxx occurs contact the service engineer and inform him of the problem and the number xxxx.

No.Error

- 61 Command is not recognised. Type I to list those control letters allowed.
- 62 No parameter group is stored with this number. Check list of parameter groups by typing L followed by pressing RETURN.
- 63 An integer number was expected as input and instead some other characters have been given.
- 64 Illegal characters given in the numeric data.
- 65 Numerical value given as input is outside the range allowed. See the Reference Card.
- 66 Instruction to jump to a parameter line cannot be followed because either the line number is not correct, or the mode has not been given so the number of lines is not defined.
- 67 An illegal window has been set: windows must be given in the format xxx-yyy where  $0 \leq xxx \leq yyy \leq 256$ .
- 68 Incorrect printout selection, illegal number given.
- 69 Incorrect printout selection, tabulation exceeds line width.
- 70 A number must follow the tab (>) character.
- 71 Illegal character in printout selection: text must be given inside apostrophes (').
- 72 Printout selection too long, more than 80 items.
- 73 Y(ES) or N(O) must be given.
- 74 Difference in position number in coding must be within the range 1 - 99. If the earlier position number is greater than the second then 1000 is added to the second and the test is made again. E.g. the difference between samples coded with position 997 and 003 is  $1003 - 997 = 6$ .

No.Error

- 75 Parameter group cannot be cleared because the group number given the second time (in response to the group number confirmation question) differs from that given first. Start clear procedure again from the READY state by typing C.
  - 80 Mode must be given when creating a new group.
  - 81 Memory size exceeded. Jump to the READY state and cancel some other parameter group from the memory.
  - 82 Cannot resume because no assay interrupted.
  - 83 Cannot resume because the mode in the parameter group has been changed.
  - 84 A checksum test has revealed an error in a parameter group in the memory. The numbers of the groups which have been lost are listed.
  - 91 One or more motors have been on too long (more than 9 minutes) so all motors have been stopped. Possible causes are:
    - a) An obstacle on the conveyor. Remove the obstacle and type letter O then press RETURN to drive the motor to the initial position.
    - b) The instrument has been left running in the A mode with no racks, or invalid racks, for more than 9 minutes. Continue as in a)
    - c) A motor is not working properly.
- Note: If the message RAM FAILURE xxxx occurs contact the service engineer and inform him of the problem and the number xxxx.