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Zeiss Z2 B-301

Démarrage

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Zeiss Z2 B-301

Startup

8. Turn on the X-Cite exacte lamp (on the left of the microscope)
9. Turn on the microscope power supply (supply 232, on the right of the microscope)
10. Turn on the microscope (on the microscope on the left side at the back)
11. Turn on the Prime camera (on the top of the top camera). **The computer must be started right after the camera.**
12. Turn on the computer
13. Log in Windows using your UdM credentials
14. Start up **Zen Blue** software

Shutdown

6. Log in the microscope log book
7. Get your sample from the microscope
8. **Clean the oil objectives with lens cleaner and lens paper**
9. Check on **FACES** if someone is using the microscope after you (same day):
If so:
 10. Inform the person that the microscope will be left on and that she/he will be responsible for turning it offIf not:
 5. Shutdown the microscope
 - a. Close **Zen Blue** software
 - b. Turn off the Prime camera
 - c. Turn off the microscope
 - d. Turn off the microscope power supply
 - e. Turn off the X-Cite exacte lamp
 - f. Cover the microscope
6. Copy your data on your external drive and **delete it from the local C: drive**
7. Log off Windows

Important Reminders

- Take back your samples including ones in the microscope
- Leave the microscope and the working area clean
- X-Cite lamp must be on **for at least 30 min** before being turned off **and vice-versa**

Storage management

- Files can be saved **temporarily** (during acquisition) on the local C: drive (desktop)
- **At the end of each session**, copy your data to your external drive and **delete it from the local C: drive**

Support

- Platform managers are available from Monday to Friday from 8:30 AM to 12:00 PM and from 1:30 PM to 5:00 PM
 - Phone: 514-343-6111 extension 5148
 - Office: D-333, Roger Gaudry Building
- Monique Vasseur, Microscopy Platform Manager, monique.vasseur@umontreal.ca
- Nicolas Stifani, Advanced Microscopy Coordinator, nicolas.stifani@umontreal.ca